Ovulation[®] test

Français	P 4
English	P 8
Espanioles	P 12
Italiano.....	P 16
Português	P 20
Protocole opératoire.....	P 23
Operating procedure	

USAGE IN VITRO UNIQUEMENT

FRANÇAIS

PRINCIPE

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ENGLISH

PRINCIPLE

The Ovulation Test is a rapid semi quantitative ELISA test on membrane with a visual reading for the determination of progesterone level in plasma or serum of dogs. It allows essentially to determine the ovulation and therefore to forecast the date of insemination or covering. Ovulation Test is an ideal complementary tool to vaginal smear examination during the heat follow up.

INTEREST

■ DETERMINATION OF THE FERTILE PERIOD :

The visible signs of heat and breeding time (such as vulval swelling, vaginal discharge, or flagging) are only an approximate indication of time of ovulation. These signs can vary by more than one week. The vaginal smear examination is necessary to confirm the oestrus or pro-oestrus stages.

The Luteinizing Hormone or LH triggers the ovulation. The day when LH increases is called Day 0. The ovulation occurs 2 days after. After ovulation, the oocytes take 2-3 days to mature, and the mature ova are alive for 48 to 72 hours. Therefore, the most fertile period for the bitch is 5-6 days after the LH peak.

According to the figure p23, the progesterone levels are low before the ovulation time (between 0 to 1 ng/ml, test bright blue). When the progesterone begins to increase, the color of the test changes from bright blue to light blue, corresponding to peak LH (day 0). Then the progesterone continues to increase the following days and the test color becomes faint blue, then white.

The pregnancy is only possible with mature ova (2 days after ovulation) and with alive sperm (5 days +/- 2). 20% of dog shows ovulation outside the period between the 10th and 14th days of heat. According to C. DUMONT (1) only two associated examinations, vaginal smears and detection of progesterone level, allow to schedule insemination or covering and to reach 95% of success with an increasing of prolificity.

■ ANOTHER INTEREST FOR DOGS :

Determination of the whelping date :

When it is used in conjunction with all other available information such as drop of rectal temperature, white milk, ultrasound examination and others, level of progesterone lower than 1 ng/ml indicates a readiness for delivery within 24 hours or 48 hours.

Another interests :

- Follow up of fertility treatment
- Preventing abortion due to luteal insufficiency by controlling the effect of adapted treatment.

■ ANOTHER INTEREST FOR BOVINE, EQUINE :

- Early pregnancy status
- Efficiency of progestagens or progestative treatment

PRESENTATION

6 tests / kit:	6 unitary test cups	1 vial of SUBSTRATE A
	6 transfer pipettes	1 vial of SUBSTRATE B
	1 vial of WASH SOLUTION 1	1 vial of Substrate Mixing Bottle
	1 vial of WASH SOLUTION 2	1 vial of Conjugate
	1 colour guide	1 pack insert

STABILITY – STORAGE

Store at 2-8°C. DO NOT FREEZE

Shelf Life : 14 months from the manufacturing date written on the box

RECOMMENDATION

■ WHEN TO TEST :

Generally 1 or 2 tests are enough if the first day of the appearance of pro-estrus is known (vaginal bleeding). The vaginal cytology can be used like a guide in order to identify the approach of LH peak and then to determine the progesterone level.

For particular cases, like infertility problem or non-regular cycle, it is recommended to tests starting from the 4th day during pro-estrus and to continue until insemination by testing every 2 days.

■ SAMPLING :

·Collect the blood sample into an EDTA or Heparin coated tube or in a dry tube. Do not use tube with serum separator gel. Do not use whole blood.

·Immediately after collection, invert the tube many times to mix it. Spin the blood down with a centrifuge or let the blood clot by staying at room temperature (2 hour) and pour the serum into a clean glass tube.

Do not use hemolysed or lipemic serum.

If the sample is not immediately used, store the sample at the fridge (2-8°C). For a long

storage, the serum or plasma should be frozen.
Label the sample.

■ **PRECAUTIONS OF USE :**

Always reseal the plastic bag after removing a test cup.
Do not exchange the test cups or the reagents from one box to another

PROCEDURE (fig. p. 23)

Take the kit out from the refrigerator at least 2 hours before using a test, in order to allow the reagents to come at room temperature. It is also possible to leave the kit a room temperature along the day.

Timing the steps 1,3 and 7 is very important.

1 With a transfer pipette, **add 8 drops of sample** to the center of the cup.
Wait 2 minutes.

2 **Add 4 drops of Wash solution 1.** Wait for the liquid to drain completely into the cup.
Repeat this step 1 time.

3 **Add 3 drops of conjugate** (red cap bottle) to the center of the cup.
Wait 1 minute.

4 **Fill the cup to the top of inner line with the Wash solution 2,** for about 20 drops.
Wait for the liquid to drain completely into cup.

5 **Prepare a fresh substrate solution** using the empty mixing bottle (blue label bottle): add one volume of each substrate A and substrate B bottles by filling their droppers up to the line. Homogenize by shaking the freshly prepared substrate.

6 **Add 4 drops of the freshly prepared substrate bottle** (blue label bottle) to the center of the cup.
Discard any unused freshly prepared substrate solution after 30 minutes as it becomes rapidly unstable.

7 **Wait 9 minutes.** Read and record the result, with the help of the colour guide.

INTERPRETATION OF RESULTS

The colour result is corresponding to the level of progesterone. Use the enclosed colour guide for evaluating and recording the result.

The following table is a guide of interpretation for the follow up of the fertile period in dog female

Color	Progesterone Level	Interpretation*
C1 Bright Blue	0 to 1 ng / ml	Baseline progesterone level. Test again every 2 days until getting a light blue colour result.
C2 Light Blue	1 to 2,5 ng / ml	LH peak Test again every 2 days in case of suspected cycle troubles.
C3 Faint Blue	2,5 to 8 ng / ml	The eggs begin to ripen after ovulation. Mate or inseminate in 1 and 3 days (and 5 days with frozen semen).
C4 White	> 8 ng / ml	Mate or Inseminate immediately. Caution: may be the fertile period is ended.

* The recommended actions above are a guideline, as no test is 100% accurate in all times and under all conditions. Therefore, the veterinarian should interpret all test result in light of the patient's clinical and/or ultrasonic examination, and historical information. This kit is intended to help the practitioner in the choice of treatment. This is the veterinarian's prerogative and responsibility. Bio Veto Test can not be held responsible for the consequences of misuse or misinterpretation of the results given by the test.

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