EN
Speed cCRP™
 Veterinary diagnostic test
 For *in vitro* use only

Fluorescence immunoassay test for the quantitative measurement of canine C-Reactive Protein (cCRP) in canine serum or plasma.

■ **CLINICAL APPLICATION**
 Canine CRP (cCRP) is a major acute phase protein primarily synthesized by the liver during the course of systemic inflammation. Marked cCRP alterations are associated with the evolution of inflammation and permit early detection as well as adapted treatment follow-up.

■ **PRINCIPLE**
 Speed cCRP™ is a rapid immunochromatographic quantitative test providing an assay of the circulating cCRP by laser-induced fluorescence analysis.

■ **PRECAUTIONS AND WARNINGS**
SAMPLE
 - Use only serum or plasma (Heparin-Li) samples.
 - Do not use tubes containing serum separator gel.
 - No interference with haemolysis or hyperlipaemia has been seen.
 - If the test is not performed within an hour, transfer the serum or plasma to a dry tube.
 - Samples can be stored in a refrigerator (+2°C to +8°C) for up to 7 days.
 - For testing beyond 7 days, the serum or plasma should be stored in a freezer (-20°C).
 - Stored samples should be allowed to reach room temperature (18°C to 27°C) and be centrifuged before analysis.
 - The presence of fibrin strands in the sample may lead to erroneous assay results due to pipette clogging and reduced volume of aspirated sample. In this case, recentrifuge the sample.

■ **PROPER PIPETTE USAGE**
 - Hold the pipette, the sample and reagent tubes vertically at all times, to ensure the correct sampling volume.
 - Use a new disposable pipette tip for each test.
 - Press and hold the pipette plunger and insert the tip into the sample, being careful not to touch the walls of the tube. Release the pipette plunger slowly while keeping the end of the tip in the sample.
 - Remove the pipette from the sample tube and check for air bubbles in the tip. If air bubbles are present, discard the sample and repeat the sampling.

■ **REAGENT AND SAMPLE MIXTURE**
 - Before using the first test of each box, register the batch calibration chip at the Speed Reader™.
 - The expiry date is indicated on the box and on each test pouch.
 - Store the tests and the reagents between +2°C and +8°C.
 - The specimen and all the materials used for the test procedure should be considered as potentially infectious and they should be disposed of in accordance with local regulations.
 - Use appropriate and clean protective equipment (gloves and gown).
 - In case of skin or eye contact with the reagent, rinse immediately with clean water and ask for medical advice.

■ **RECOMMENDATIONS**
 - Before using the first test of each box, register the batch calibration chip at the Speed Reader™.
 - The expiry date is indicated on the box and on each test pouch.
 - Store the tests and the reagents between +2°C and +8°C.
 - The specimen and all the materials used for the test procedure should be considered as potentially infectious and they should be disposed of in accordance with local regulations.
 - Use appropriate and clean protective equipment (gloves and gown).
 - In case of skin or eye contact with the reagent, rinse immediately with clean water and ask for medical advice.

■ **TEST PROCEDURE**
 Speed cCRP™ cannot be visually interpreted and should be read only with a Speed Reader™ analyser. For complete instructions, please refer to the Speed Reader™ analyser operator's guide.

- Do not mix test devices, reagent tubes and calibration chips from different batches.
 - Allow the reagent and the test device to reach room temperature (18°C to 27°C) for at least 30 minutes before use.

■ **MATERIALS FOR EACH TEST**
 1 test device, 1 reagent tube, 1 dropper cap, 1 pipette tip and the 10µL pipette

■ **SPEED READER™ PREPARATION**
 - Press "New Test" on the main screen before starting the test procedure.
 - Follow the on-screen instructions.

■ **SAMPLE TRANSFER**
 - Using the tip of the dropper cap, perforate the aluminium foil in order to open the reagent tube.
 - Place a new disposable tip on the pipette.
 - Transfer 10 µL of the sample into the reagent tube, following the instructions for correct pipette usage.
 - Carefully seal the reagent tube with the dropper cap.
 - Gently mix the contents by inverting the tube at least 5 times. DO NOT SHAKE.

■ **SAMPLE APPLICATION**
 - Remove the plastic cover from the dropper cap and discard the first 2 drops of the mixture, to wash out the dead space of the dropper cap.
 - Completely remove the cartridge from Speed Reader™ and add 2 drops of the mixture into the sample well of the test device.

■ **READING**
 - Insert the test device into the Speed Reader™'s cartridge holder, when prompted by the analyser, in the direction indicated by the arrow (the sample well first) and follow the on-screen instructions.
 - Read the result on the screen.

■ **INTERPRETATION**
 cCRP concentration is measured in mg/L, between 5 mg/L and 250 mg/L

Reference values:	
Normal: < 15 mg/L	=> Systemic inflammation is unlikely
Borderline: 15-40 mg/L	=> Mild systemic inflammation => Progressive systemic inflammation (onset or in course of resolution)
High: > 40 mg/L	=> Consistent with systemic inflammation

These recommendations are for guidance only. The aim of this test is to measure canine CRP blood concentrations. Interpretation of the result by the veterinarian should always take into account the history, clinical examination, and any further diagnostic test results, as no diagnostic method is 100% accurate. The definitive diagnosis is the prerogative and responsibility of the veterinarian. Bio Veto Test and its distributors cannot be held responsible for any consequences linked to incorrect use of this test or misinterpretation of the results.

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