# WITNESS™ | Giardia

## CANINE AND FELINE GIARDIA ANTIGEN TEST KIT.

### WITNESS™ | Giardia

ENGLISH

#### I. GENERAL INFORMATION

Giardia intestinalis (syn. G. duodenalis or G. lamblia), is a gastrointestinal protozoan parasite of various species including dogs and cats. Giardia can cause a variable clinical picture ranging from asymptomatic carrier to acute or chronic diarrhea. This flagellated parasite has two different morphologic stages: an axenial intestinal multiplication form (trophozoite), and an encysted form, that is excysted and is readily infectious (cyst). Infection occurs via oral-fecal contact with an infected animal or ingestion of contaminated water or material from the environment. Symptoms appear around 3 weeks post-infection, with a normal incubation period of around 8 days. Typical clinical manifestations include diarrhea, weight loss, and growth retardation. The disease may be self-limiting or lead to chronic diarrhea and intestinal malabsorption. This test can be used as a screening tool in young animals such as dogs and cats aged from between 6 weeks and 5 months. It can also be used in older animals suffering from chronic diarrhea. Giardia asymptomatic carriers may exist, test results and consequent treatment decisions should always be considered in the context of all available clinical information.

#### II. TEST PRINCIPLES

This test uses red polystyrene micro-spheres coated with monoclonal antibodies to *Giardia intestinalis* antigen. Blue microspheres are also present as a test control. Positive antigen present in the faecal sample reacts with the latex particles coated with monoclonal antibodies specific to the antigen. This particle complex of latex/antibodies/parasite antigen subsequently migrates through the test reaction zone via chromatography. In this zone, anti-Giardia antibodies react with the particle complex of latex/antibodies/parasite antigen resulting in a red/pink line appearing (positive test).

#### III. SAMPLE COLLECTION

- Test material required: canine faeces, feline faeces

#### IV. SAMPLE STORAGE

Faecal samples may be stored at +4°C for 24 hours. If longer storage is required, the samples should be frozen at -20°C. Faecal samples must be allowed to reach room temperature before testing.

#### V. KIT CONTENTS

- 5 pouches, each containing 1 test device and desiccant.
- 5 disposable sample extraction pipettes.
- Extraction buffer bottle.
- Instructions for use.


#### VI. PRECAUTIONS

- Do not use components after expiration date on the pouch.
- Store-test kits between +2°C and +25°C. Do not freeze.
- Use the test within 10 minutes of opening the sealed pouch.
- Avoid touching or damaging the membrane at windows (1, 2) or (3).
- The device should be placed on a flat, horizontal surface for the duration of the test.
- Use a separate extraction pipette for each sample.
- Hold extraction pipette vertically when dispensing sample into the test well.
- Handle all reagents and samples as biohazardous material.
- For veterinary use only.

#### VII. TEST PROCEDURE


##### 1. SAMPLE COLLECTION

- Remove the top, containing the cotton swab, from supplied sample extraction pipette.
- Dispense Extraction Buffer from the bottle supplied into the base of sample extraction pipette until the buffer reaches the graduation mark.
- Holding the top of the sample extraction pipette, coat the cotton swab with faeces from a stool sample, or from an anal swab. It is critical that the sample contains between 60mg to 100mg of faeces; so ensure the cotton part of the swab is entirely covered with faeces.
- Insert the swab into the dispersed extraction buffer in the base of the sample extraction pipette. Securely fit the top onto the base, then tap the pipette to mix the contents for 5 seconds ensuring good sample extraction.

##### 2. SAMPLE TESTING


- First, tear open the pouch provided and place the test device on a flat horizontal surface for the duration of the test, to allow correct migration.
- Then snap the upper part of the sample extraction pipette at the blue line. This contains the extracted faecal sample.
- Invert and hold the sample extraction pipette in a vertical (upright) position and gently squeeze the base to dispense five drops of sample onto the sample well (window (1)). Avoid touching the well itself.
- If fluid is not seen to migrate past window (2) within one minute, add one additional drop of sample into the sample well, window (1).

##### 3. READING TEST

- Wait 5 minutes, and then look for the presence or absence of a red/pink band in window (2) and a blue band in window (3).
- Sample test results are read in window (2). The control test band is read in window (3).
- Note:**

  - A control blue band may appear in window (3) before 5 minutes has elapsed, however the test may not have finished. A red/pink band in window (2) can develop more slowly than the blue control band in window (3) so 5 minutes should elapse before reading.
  - The test may be completed and read before 5 minutes only if bands are visible in both windows (2) and (3).

##### VIII. RESULTS

- Negative:** A sample is negative for *Giardia intestinalis* antigens when there is no red/pink band visible in window (2) whilst a blue control band is visible in window (3).
- Positive:** A sample is positive for *Giardia intestinalis* antigens if there is a red/pink band visible in window (2) and a blue control band is visible in window (3).
- Note:**

  - A test is only valid if a blue control band is present in window (3). If there is no blue control band visible in window (3), the test is invalid and should be repeated.

##### Remember:

- A test result should always be interpreted in the context of all available clinical information and history for the dog or cat being tested.

##### TECHNICAL ASSISTANCE

If you need assistance, please contact your local representative.

##### SYMBOL DESCRIPTIONS

	Use by (expiration date)		Consult instructions for use
	Lot		In vitro diagnostic medical device
	SN		Manufacturer
	Temperature limitations (storage temperature range)		

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- Négatif :** A sample is negative for *Giardia intestinalis* antigens when there is no red/pink band visible in window (2) whilst a blue control band is visible in window (3).
- Positive :** A sample is positive for *Giardia intestinalis* antigens if there is a red/pink band visible in window (2) and a blue control band is visible in window (3).

##### Note :

- A test is only valid if a blue control band is present in window (3). If there is no blue control band visible in window (3), the test is invalid and should be repeated.

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