BOVINE VIRUS DIARRHEA ANTIGEN TEST KIT

WITNESS™

BVDV

ENGLISH

I. GENERAL INFORMATION

WITNESS® BVDV detects Bovine Viral Diarrhea Virus (BVDV) Erns antigen in bovine ear tissue samples.

II. TEST PRINCIPLES

WITNESS® BVDV is a simple test based on Rapid Immuno Migration (RIM®) technology. BVDV Erns antigen is detected in ear notch samples from persistently infected (PI) animals. Sensitized particles bound to the BVDV Erns antigen present in a sample from BVDV infected animals migrate along a membrane. The complex is then captured on a sensitized reaction line where its accumulation causes the formation of a clearly visible pink/red band. A control band on the strip ensures that the test was performed correctly.

III. SAMPLE COLLECTION

Before collecting the ear notch, dip the notching tool in disinfectant and then rinse away residual disinfectant with copious quantities of clean water. Select a clean portion of the animal's ear and avoid blood vessels to harvest a sample. Refer to the table below to obtain a proper ear notch size.

	Large ear notch (≥ 1 cm)	Small ear notch (≈2-3 mm and < 1 cm)
SIZE REQUIRE- MENT	Ear notch is at least 1 cm on one side	Ear notch is at least 2-3 mm in diameter and <1 cm in all dimensions

IV. SAMPLE STORAGE

For short term storage, ear notches should be stored refrigerated at 2 - 7 °C.
If refrigerated storage is not available, ear notches can be stored (up to 3 days at 15 - 30 °C).
For long term storage, the ear notch samples should be kept frozen (-20 °C or colder).
Allow samples to come to 15 - 30 °C before testing.

V. KIT CONTENTS

- 1 Desiccant vial containing 10 test strips
- 2 Bottles Extraction buffer (10 mL each)
- 10 graduated tubes with caps
- 10 Test tubes
- 10 Disposable pipettes
- 1 Rack with holders for 5 test tubes

Instructions for use

VI. PRECAUTIONS

- Do not use this kit or any of its components after the expiration date.
- Kit should be stored at 2 - 25 °C. Kit should not be frozen.
- The vial holding the test strips contains a desiccant and should be kept tightly closed when not in use.
- Use the test immediately after removal from the desiccant vial (within 10 minutes).
- Test strips should only be handled in the upper, labeled region. Avoid contact with the surface of the test strip.
- Use a separate pipette for each sample.
- Hold the pipette vertically when dispensing buffer from the extracted sample.
- Extraction buffer is preserved with sodium azide.
- For animal use only.


VII. TEST PROCEDURE

1. SAMPLE PREPARATION


	Large ear notch (≥ 1 cm)	Small ear notch (≈2-3 mm and < 1 cm)
EXTRACTION BUFFER	1.5 mL	0.5 mL
MIX BY INVERSION	4 times	4 times
INCUBATION TIME (at +15 - 30 °C)	At least 10 minutes (10 minutes to 1 hour)	At least 1 hour (1 hour to 3 days)
EXTRACTED SAMPLE STORAGE	If the extracted ear notch samples will not be tested immediately they should be stored at 2 - 7 °C for up to 48 hours. For prolonged storage, samples should be kept frozen (< -20 °C).	

2. BUFFER DISPENSING

- First add the appropriate volume of extraction buffer to the graduated tube as follows:
 - 1.5 mL of extraction buffer for a large (≥ 1 cm) ear notch
 - 0.5 mL of extraction buffer for a small (≈ 2-3 mm and < 1 cm) ear notch


- Invert the tube containing the ear notch and buffer 4 times to mix, make sure the ear notch is immersed in the buffer.
- Incubate the ear notch samples within the extraction buffer at 15 - 30 °C for the appropriate length of time:
 - Incubate the large ear notch samples for at least 10 minutes
 - Incubate the small ear notch samples for at least 1 hour

3. TESTING SAMPLES

- Use a disposable pipette to dispense 250 µL of the extracted sample into the test tube provided. (Note: To provide approximately 250 µL of extracted sample, insert the pipette into the extraction solution, squeeze the bulb, draw up the liquid from the extraction tube until the liquid reaches the nearest graduated line and then dispense all the liquid into the test tube.)
- Grasp only the labeled portion of the test strips to remove one strip from the desiccant vial for each sample to be tested.
- Place the test strip directly into the test tube so that the red pad on the test strip is submerged in the extracted sample.

- Incubate the test strip in the extracted sample for 15 minutes.
- Remove the test strip from the test tube.
- Read results immediately.

VIII. RESULTS

Reading Results

- After 15 minutes, observe the presence or absence of pink/red bands in the center of the test strip between the two absorption pads.
- The control band appears in the upper end of the test strip (closest to the labeled region), while the sample test results are read in the lower part of the test strip.
- Discard the test strip and test tubes in an appropriate biohazard container in accordance with local regulations.

Valid Results

The test is **VALID** if a pink/red band develops in the upper Control part of the test strip. The absence of the Control band indicates that the test is invalid and must be repeated.

Interpretation of Results

- **POSITIVE** for BVDV: Two pink/red bands (Control band and Test band) are clearly visible on the test strip (C & T). With WITNESS® BVDV, as with other testing methods, a positive result may not always indicate the animal is persistently infected; the animal might be acutely infected with BVDV. The same animal should be retested after three weeks to confirm its PI status. If the result remains positive, the animal will be identified as PI animal.
- **NEGATIVE** for BVDV: A single pink/red band is present in the upper part of the test strip (C).

Note: The test is complete and may be read before 15 minutes if both control and test pink/red bands are visible. The presence of only the control band on the strip before 15 minutes does not mean that the test is complete, as a test band may appear more slowly than the control band.

SYMBOL DESCRIPTIONS

	Use by Date (expiration date)		Authorized Representative in the European Community
	Batch Code		Consult Instructions for Use
	Serial Number		In Vitro diagnostic medical device
	Temperature limitations (storage temperature range)		Manufacturer

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