Dermatophytest™

English	1
Français	4
Español	7
Português	10
Italiano	13
Ελληνικά	16
Deutsch	19
Netherlands	22

FOR *IN VITRO* USE ONLY

ENGLISH

■ CLINICAL APPLICATION

Dermatophytosis is a fungal infection of the skin encountered in dogs, cats, rodents, and lagomorphs. It is essentially caused by various genera of dermatophytes (fungi) and notably *Microsporum* and *Trichophyton*. Microscopic examination of hair and skin scrapings and examination under a Wood's lamp do not provide a definitive diagnosis. Only the growth of dermatophytes from a sample of suspected material constitutes a final diagnosis.

Dermatophytosis is a zoonosis.

■ PRINCIPLE

Dermatophyttest is a selective agar medium for the diagnosis of Dermatophytosis. When pathogenic dermatophytes of the genus *Microsporum* and *Trichophyton* are grown at room temperature, they have the particularity to metabolize nitrogen compounds of the medium. Thus, dermatophytes produce nitrogenous wastes which alkalize the agar medium and lead to a color change of the agar from yellow-orange to red in less than a week.

■ PROCEDURE

1/ PREPARATION OF THE LESION:

- Wash the skin lesion with 70° alcohol or soap (non-fungicidal).
- Dry with absorbent paper.
- This eliminates any saprophytic agents that could cover the growth of dermatophytes on the agar medium.

2/ SAMPLE:

- **Skin scrape:** scrape the lesion with a glass slide or the blunt side of a scalpel blade, from the edge of the infected area towards normal skin. Unload hairs and scales onto the agar medium.
- **Hair pluck:** deposit hairs and scales onto the agar medium.
- **Swab:** moisten the extremity of a sterile swab with some drops of moistening solution. Wipe the swab over the entire lesion and then discharge the swab onto the agar medium using back and forth movements, without crushing the agar medium.

3/ CULTURE:

- Record the name of the animal and the date on the dish containing the agar medium.
- **Do not close the agar dish fully** during incubation; culture is facilitated by an exchange of air with the outside.
- Allow to grow at room temperature, away from the light.

4/ READING:

- Read two days after inoculation and then every 2 days until 14 days.

- **The growth of dermatophytes induces a colour change in the agar from yellow-orange to red.**

- **POSITIVE TEST:** slight red colour change associated with the presence of clear, downy, or filamentous colonies.

A colour change due to the growth of saprophytic microorganisms is slower to appear, appearing 1 to 2 weeks after that induced by dermatophytes. Furthermore, the morphology of these saprophytic colonies is not characteristic (deeply pigmented filamentous colonies or mucous colonies).

■ RECOMMENDATIONS

• SHELF LIFE:

- 16 months between +2°C and +8°C from the manufacturing date.
- The expiry date is printed on each box. Keep the test away from sources of excessive heat or cold.
- If the agar medium is stored between +2°C and +8°C, it does not need to be brought to room temperature before inoculation.

• PRECAUTIONS:

- Protective gloves should be worn when handling samples and inoculating the agar.
- Avoid incubation in the light; dermatophytes prefer to grow in a dark environment.
- Do not close the agar dish fully during incubation.
- After inoculation, discard contaminated instruments or clean with alcohol before sterilisation.
- All used agar media must be correctly sealed and disposed of in the hazardous waste bin to prevent the dispersion of dermatophytes.

• OTHER RECOMMENDATIONS:

The above recommendations are only guidelines; no test is 100% accurate in all times and under all conditions. The purpose of this test is to assist the practitioner with the diagnosis of dermatophytosis. All test results should be interpreted in the light of the patient's history, clinical examination, and the results of any other further diagnostic tests. The definitive diagnosis remains the veterinarian's prerogative and responsibility.

Bio Vétro Test and its distributors cannot be held responsible for the consequences of misuse or misinterpretation of the results given by the test.

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