Speed[®] BIOGRAM

Françaisp. 3

Englishp. 10

USAGE IN VITRO UNIQUEMENT

BREVET Bio Vétro Test

INTERET CLINIQUE

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Bio Veto Test Patent

CLINICAL INTEREST

Speed® Biogram is a diagnostic tool for the identification of pathogenic bacteria (and yeasts) responsible of dogs and cats dermatological and bacterial (and fungus) infections, and urine or auricular infections.

It allows obtaining a quick antibiotic susceptibility testing among the veterinarian therapeutic antibiotics available.

This test, realizable in a veterinarian structure, is obtained directly from a liquid or a cellular sample. Within 24 hours, it suggests you the specific antibiotics to use for the pathogen bacterial agent, and identifies the bacteria within 48 h (fungus in several days).

This rapid patented technology takes into account the synergistic or antagonist effects of the different pathogen agents, of the infected environment factors (direct antibiogram) and the concentration of the infected site germs (inoculum effect)

PRINCIPLE

Speed® BIOGRAM is a micro-tray culture containing:

- 13 wells for sensitivities and bacterial resistance for 13 veterinarian antibiotic molecules.
- 6 wells for the identification of the pathogen(s) bacteria.
- 2 wells for the identification of yeasts (malassezia / candida)
- 3 wells (BVT Patent) which points out the reading delay to consider the pathogen agents' antibiogram (and not the antibiogram of an eventual non pathogen agent)

One of these 3 wells is called I.L well (Incubation Limit), it is the only well that doesn't receive the sample. This control well is considered as a signal, which limits the reading time of the antibiogram. I.L well reading is made concurrently to a second +/- well, which receive the sample. +/- Well is a reference solution for the growth of the germs, which color change corresponds to a bacterial concentration 10^3 UFC/ml, that is to say an infection.

Results obtained with Speed® BIOGRAM are read well by well. Bacterial growth is revealed by a color change. It shows you a resistance of the pathogen agent to the corresponding well antibiotics (concerning the 13 antibiotic wells), and indi-

cates the type of pathogen agent present (concerning the 8 bacterial identification wells).

PACKAGING

5 tests kit ref. SB 5 T

- 5 packed Speed® BIOGRAM trays with individual adhesive labels.
- 5 vial preservative mediums
- 5 growth medium vials
- 1 "Supp Staph" vial
- 1 paraffin oil vial
- 5 sterile pipettes
- 5 reading carriers for BVTuve (incubator)
- 5 color trays for customers
- 1 operating procedure

1 BVTuve warning plate (ref: BVU) to incubate Speed® BIOGRAM at 37 ° can be provided upon request.

OPERATING PROCEDURE

■ SAMPLING

Every cutaneous, urinary or auricular sample can be used, only if it is obtained from sanitized environment and if it doesn't contain antiseptics (Iodine, Alcohol, Phenol...) or local antibiotics.

Discharge immediately and delicately the sample in the kit's preservative medium. Homogenize the medium.

You can store this vial at ambient temperature (or keep it 48h between 4 and 8 °C) until the deferred use of the Speed® BIOGRAM kit.

■ PROCESSING

1 - With a new tray just taken out from its pouch, identify the animal with its name and date, then take off the adhesive strip and paste it on one side of the tray in order to conserve the wells identification.

2 - Take a medium culture vial, add 3 drops of this vial in the I.L well with a pipette (which well does not receive the sample)


3 - Discharge the cell sample or 2 drops of biological liquid with the pipette in the preservative medium, and then homogenize the mix.

4 - Transfer 4 drops from the preservative medium to the opened medium culture and homogenize

- 5 - Inoculate the tray's wells one by one, adding with the same pipette 3 drops of this medium culture, mixed in each well, except the I.L well.
- 6 - Add in the Staph well 2 drops of the " Supp. Staph. ".
- 7 - Add 2 drops of paraffin oil in each well, except in E. coli and pseudo wells
- 8 - Stick adhesive label back on the tray, and put the tray on its reading carrier. Put the reading carrier at the center of the plugged BVTuve, in a transversal sense in order to harmonize temperature conditions in the whole tray. (This sense is better than the longitudinal one, parallel to both sides of the BVTuve.)
- 9 - Use BVTuve in a room where the temperature is controlled and homogeneous. (The ideal is between 20 and 25°C)

■ **READING**

1 - Antibiogram reading: begin the reading of the antibiotic profile 18 to 24 hours later, when the color of the I.L well is turning from red to yellow. The I.L. well turns to yellow and the growth reference +/- well turns from colorless to red. This color change of the +/- well is typical of an infection, for a pathogen bacteria concentration \geq to 10^3 UFC/ml. In this case, an antibiogram profile reading can be made.

Identify sensitivities and resistances well by well and for each antibiotic as described:




Table of Speed® BIOGRAM antibiotics

Antibiotic molecules	
Amoxicillin	Clindamycin
Amoxicillin + clavulanic Ac.	Flumequin
Cefalexin	Spiramycin
Enrofloxacin	Gentamicin
Marbofloxacin	Neomycin
Difloxacin	Polymixine B
Doxycyclin	

2 - Reading of the bacterial fungus identification: read the wells individually after 24 to 48h, in order to identify a color change or a bacterial growth. Only the *Malassezia* well can be read after 3 to 5 days of incubation.


Streptococcus: surface blue color change


Pseudomonas: blue color change


Staphylococcus: yellow color change
(read at least after 48 hours)


Proteus: fushia red color change
+ Color change of the enterobacteria well


Enterobacteria: yellow color change


Candida: Yellow color change

Escherichia coli: surface blue color change
+ Enterobacteria yellow change of color

Malassezia: Black color change
(Between 3 to 5 days)

Reading associations are possible.

Candida and *Malassezia* yeasts do not show any antibiotic profile on the tray.

■ INTERPRETATION

For each reading of susceptibility and resistance profile of the searched bacteria, a sheet of results allows you to register them in order to transmit the whole antibiogram to the animal's owner, associated with your treatment choice. Antibiotics are adapted to precise dermatological, urinal or auricular indications, which you must choose according to antibiogram results and to the type of the sample.

RECOMMENDATION

During the well reading, we recommend to let the tray on the BVTuve. If the reference growth +/- well does not turn to red after 18 to 24 h of incubation, let incubate it more during 18 to 24h, in order to confirm the lack of pathogen germs in the sample.

SHELF LIFE - STORAGE

14 - months shelf life at 4-8°C from the manufacturing date written on the kit box. Before using, we recommend to let the whole reagents and the tray at least 15 min. to room temperature, and to check the functioning of the BVTuve. We recommend you not to use reagents vials from various batches.

All batches are QC controlled and satisfied the necessary requirements concerning in-vitro diagnostics.

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